# Monitoreo bacterial

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# Agenda

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1. While in the field, correctly label 1 Whirl-pak® bag with a permanent marker for the blank/control.
2. Put on latex gloves and remove the perforated seal from the top of the Whirl-pak® bag.  
**IMPORTANT!** Do not touch the inside of the Whirl-pak® as this will contaminate your sample and alter the results.
3. Use the two small white tabs to pull open the bag.
4. Fill the Whirl-pak® bag 2/3 full with distilled water.
5. Grab the ends of the twist ties and “whirl” or spin the bag tight. Cross the twist ties to close the bag.
6. Make sure the bag is closed securely by inverting the bag to test the seal (no water leaks out).
7. Immediately place the Whirl-pak® bag into a properly disinfected cooler with ice and store there throughout your sampling event.

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1. Correctly label new Whirl-pak® bag with a permanent marker for the sample/site information.
2. Put on latex gloves and remove the perforated seal from the top of the Whirl-pak® bag.  
**IMPORTANT!** Do not touch the inside of the Whirl-pak® as this will contaminate your sample and could alter the results.
3. Use the two small white tabs to open the bag.
4. While holding the yellow twist ties place the bag in the water at mid-stream, mid-depth or in a well-mixed area and allow the water to flow into the bag. Fill the bag with water up to 2/3 full.  
\*Remember to collect the water sample (at least wrist deep) upstream of where you are standing. If sampling from a bridge, use a rope tied to a small disinfected bucket to grab the sample.
5. Grab the ends of the twist ties and “whirl” or spin the bag until tight. Cross the twist ties to close the bag.
6. Make sure the bag is closed securely by inverting the bag to test the seal (no water leaks out).
7. Immediately place the Whirl-pak® bag into a cooler with ice.
8. Optimal holding time for samples on ice or refrigeration is less than 6 hours but no more than 24 hours.
9. Properly dispose of gloves.

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Photo credit: By Sandypichapa - Own work, CC BY-SA 4.0,  
<https://commons.wikimedia.org/w/index.php?curid=49356166>

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1. Clean working area with disinfectant spray and let dry.
2. Put on latex gloves. NOTE\* You should always wear these when handling the plates (even when going to read them).
3. Correctly label plates (1 for the blank & 3 for the site sample), and lay them on a clean, flat surface. Plates should indicate stream name, site number and the incubation start time and date. See below figure for examples of how to label plates.


Photo credit: Georgia Adopt-A-Stream directions for Bacterial Monitoring

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4. Gently shake Whirl-pak® bag to ensure an even mix of sample.
5. Place the Whirl-pak® bag in a cup to keep from spilling and open the bag using the white tabs.  
**IMPORTANT!** Do not touch the inside of the Whirl-pak® as this will contaminate your sample and could alter the results.
6. Carefully remove pipette tip from sterile container. Don't touch the pipette tip inside of the sterile container and practice caution to ensure that the tip is not contaminated thereafter.
6. Pipette 1 ml of the sample using the fixed-volume pipettor.
7. Lift the top film of the Petrifilm™ plate and dispense 1 ml of sample on the center of the circular plate.
8. Slowly roll the top film down onto the sample until the plate is completely covered to prevent trapping air bubbles. Do not touch the center of the petrifilm plate.
9. If necessary, distribute the sample evenly by using the 3M® spreader or slightly tilting the Petrifilm™ plate back and forth. Tilting too much will cause the sample to pour out of the plate.
10. Leave plate undisturbed for one minute to allow the gel to solidify and then place in the incubator.
11. Repeat: Plate two more samples for a total of three plates per sample site.

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Photo credit: By Uncle Zooey - Own work, CC BY-SA 4.0,  
<https://commons.wikimedia.org/w/index.php?curid=48998391>

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1. Plan to turn on the incubator prior to plating to ensure it will be ready. Place the incubator lid on top.
2. Insert the thermometer into the incubator.
3. Once the incubator is at  $35^{\circ}\text{C} \pm 1^{\circ}\text{C}$ , place the processed Petrifilm plates in the incubator and reset the thermometer.
4. Incubate plates in a horizontal position, with the top film side up, in stacks of up to 20 plates. **Incubate for  $24 \pm 1$  hour at  $35 \pm 1$  degrees Celsius.**
5. After 24 hours, remove plates (with gloves on) and count *E.coli* colonies.
6. Record the minimum and maximum temperatures that are displayed on the thermometer after incubating, as well as the time in/out of the incubator.
7. Record all data on the Bacterial Data Form.
8. Dispose of plates by spraying them with an appropriate disinfectant and placing in a sealed zip lock bag, discarding in the trash.

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